

AESTHETE MEDTECH

FocusShock

Focused shock wave therapy

15-inch	1–20 Hz	7 silicone pads
Capacitive touchscreen	Pulse frequency	Replaceable interfaces

A compact therapy platform.

FocusShock combines a touchscreen console with a connected treatment handpiece and replaceable silicone pads. Controls on the handpiece provide access to frequency, intensity and session functions.



Specifications

Display
15-inch capacitive touchscreen

Pulse frequency
1–20 Hz

Intensity setting
1–10 bar

Power
10–500 W

Operating voltage
AC 100–240 V, 50–60 Hz

Main unit dimensions
39 × 43 × 43 cm

Transport case
50 × 56 × 66 cm

Gross weight
26.2 kg

Integrated organisation

Side-mounted holders keep the handpiece and silicone pads close to the console.

Touchscreen console / Connected handpiece / Replaceable silicone pads

Control at the handpiece.

Frequency and intensity can be adjusted from the handpiece interface. A local display, menu keys and start/pause control bring key functions into the operator’s hand.



Frequency

Adjustable pulse frequency across a 1–20 Hz range.

Intensity

A separate strength setting adjusts the output level.

Session control

Menu access, start/pause and an operation/stop button.

Seven silicone pads

A choice of replaceable pad profiles supports different contact configurations.



Pad 1 **Pad 2** **Pad 3** **Pad 4** **Pad 5** **Pad 6** **Pad 7**

03 / APPLICATIONS

Targeted treatment.
A versatile platform.

Treatment Areas

A handpiece-based format for localised musculoskeletal treatment.

Shoulder & elbow

Upper-limb treatment areas.

Back & hip

Trunk and hip treatment areas.

Knee & lower leg

Lower-limb treatment areas.

Ankle & foot

Distal lower-limb treatment areas.

Purposeful design

A large touchscreen, handpiece controls and organised accessory holders keep the platform compact and accessible.

AESTHETE GH LLC

30 N GOULD ST STE R

SHERIDAN, WY 82801, USA

Email: info@aesthetegh.com

Applications



Shoulders

Relieves shoulder pain and improves mobility.



Back

Reduces lower back pain and promotes recovery.



Knee

Relieves knee pain and enhances joint function.




Elbow

Reduces inflammation and accelerates recovery.



Ankle

Relieves ankle pain and supports faster healing.



Legs

Relieves muscle pain and enhances performance.



Hip

Relieves hip pain and improves flexibility.



ED function treatment

Promotes blood circulation and improves sexual function naturally.

Treatment Parameters

Treatment area	Silicone pads	Impacts	Frequency (Hz)	Intensity (bar)
Shoulder Calcium area of aponeurosis	4, 5, 6	1500–2000	4–5	4–7
Scapulohumeral periarthritis	4, 5, 6	1500–2000	4	3–5
Elbow Tennis elbow / Golfer’s elbow	2, 3	1500–2000	4–5	2–5
Wrist Carpal tunnel syndrome / tenosynovitis / TFCC syndrome	1, 2	1500–2000	4	2–5
Thigh Iliotibial band syndrome / hamstring strain	5, 6, 7	1500–2000	4	2–5
Knee Patellar chondromalacia / patellar tendinitis	2, 3	1500–2000	4	2–5
Front calf Tibial anterior muscle syndrome / tibial medial stress syndrome	1, 2	1500–2000	4	2–5
Rear calf Gastrocnemius muscle spasm	3, 4, 5	1500–2000	4	2–5
Foot Achilles tendinitis	2, 3	1500–2000	3	2–6
Foot Calcaneal bone spur	3, 4	1500–2000	4–5	4–7
Foot sole Plantar fasciitis	2, 3, 4	1500–2000	4	3–7